



P H O E N I X 4K

ULTRA HIGH DEFINITION
MINIMALLY INVASIVE SURGERY

4K
ULTRA HD


 rutledge



P H O E N I X
 4K ULTRA HIGH DEFINITION
 MINIMALLY INVASIVE SURGERY



Embark On A New Era In Animal Medicine With Our Advanced Endoscope Featuring 4K High-Definition Imaging. Immerse Yourself In The Unparalleled Clarity And Precision Of This Technology, Offering Veterinarians An Incredibly Sharp View Of Internal Structures During Diagnostic Procedures. The Superior Resolution Of 4K Not Only Enhances Your Ability To Detect Subtle Abnormalities But Also Provides A More Detailed And Comprehensive Understanding Of The Animal's Health, Revolutionizing Your Approach To Informed Decision-Making.



Incorporating 4K High-Definition Into Your Veterinary Practice Isn't Just About Seeing More – It's About Seeing Better. Elevate Your Diagnostic Capabilities And Redefine Your Approach To Animal Care With Our State-Of-The-Art Endoscope. Delivering Unparalleled Image Quality, Our 4K Technology Sets A New Standard In Animal Medicine, Empowering Veterinarians To Provide The Highest Level Of Care For Their Patients.

System Includes

- PHOENIX 4K UHD Camera System**
- PHOENIX 4K LED 50W LIGHT SOURCE**
- PHOENIX 4K 40L Co2 Insufflator**
- PHOENIX 4K LAPAROSCOPE**
- PHOENIX 4K ARTHROSCOPE/RHINOSCOPE**
- LG Healthcare 55" 4K Surgical Display**



The Phoenix 4K UHD System Includes 3 Degrees Of Image Deflection for Each Of The Endoscopes, 0, 30 and 70 Degree Angles. Rod Lenses That Are Finely Sculpted And Polished For 4K Resolution.

Using The Phoenix 4K Endoscopy System Will Put Your Practice On The Same Level Of Care As Any Human Hospital or Orthopedic Practice.

For Better Patient Outcomes.



Endoscope Model	Diameter	Total Length
PS-550L A/B/C Laparoscope	5.5mm	308mm
PS-100L A/B/C Laparoscope	10.0mm	308mm
PS-400A/B/C Arthroscope	4.0mm	157mm

Peripheral Unit Model	Primary Feature	Secondary Feature	Universal Features
SC-4K Camera Processor	4K Wide-Gamut Color	USB-A	C-Coupler
LS-4K LED Light Source	50W Dual LED	Bulb Redundancy	Olympus Style Light Guide
SI-4K Co2 Insufflator	50L Per Minute Flow Rate	Digital Pressure Control	Touch Screen Controls





P H O E N I X 4K
ULTRA HIGH DEFINITION
MINIMALLY INVASIVE SURGERY

A Full Line Of Surgical Tools And Accessories...



To Order or For More Information:

Call Us @ +1 (888) **ENDOVET**

Online @ endoscopes.vet

Socials @[thescopeguys](https://www.instagram.com/thescopeguys)

Or Call Your Rutledge Medical Representative

