




AURORA XP



rutledge



**Ben Buchannan, DVM, DACVIM - using the Aurora XP-1P First Generation*



1080p High-Definition Endoscopy

Total Endoscopy Solution For Veterinary Medicine...




A Complete Transformation For Animal Health With The All New AURORA XP Endoscopy System. This Leap Forward In Technology Allows Practitioners To Use One Tower For All Procedures Using A Scope. Both Flexible and Surgical/Rigid. Plus, For Vets Who Need Mobile Capabilities, We Offer A Processor Just For You! All AURORA XP Systems Come Equipped With Video And Image Capture Built-In To The Video Processor Units. Capture Those Images Quickly And Conveniently From The Control Body Of The Flexible Scope Or Camera Head. No Distractions.



AURORA XP-2P

Fully Portable Endoscopy Is Now Possible! And With Our Upgrades From The XP-1P, This Will Become Your Workhorse. A 40W LED Light Source Backed By 1080p HD Video Imaging With Improved Image Processing. And Now Patient Data Can Be Displayed With The XP Keyboard.



AURORA XP-2


Our 50W Dual LED Light Source Coupled With Additional Features For In-Clinic Procedures. Compatible With All AURORA XP Video Endoscopes.





1080p High-Definition Endoscopy

XP Series Flexible Videoscopes



Endoscope Model	Outer Diameter	Instrument Channel	Working Length
XP-VS3300	12.9mm	3.7mm	3.3M

Endoscope Model	Outer Diameter	Instrument Channel	Working Length
XP-VS1500	8.3mm	2.8mm	1.5M

Endoscope Model	Outer Diameter	Instrument Channel	Working Length
XP-VS1200	6.3mm	2.0mm	1M




1080p High-Definition Endoscopy

XP Series Rigid Videoscopes



1080P
FULLHD



Also Works As...

PERIOCAM


High-Definition Equine Dentistry

Endoscope Model	Outer Diameter	Angle Of View	Length
XP-VS5500R	5.5mm	0°	330cm
XP-VS550T	5.5mm	30°	330cm



1080p High-Definition Endoscopy

XP Series Accessories





**Hadley Hill Equine - Detroit, MI - using the XP-1P First Generation*



AURORA XP



2828 E 51st St
Suite 302
Tulsa, Oklahoma 74105

endoscopes.vet
+1 (888) **ENDOVET**
sales@endoscopes.vet
@thescopeguys - FB, X, IG