


**Get The  
Diagnosis  
This Fast.**



**AURORA**  
Xpress



rutledge



## Built For The Track...

Introducing our cutting-edge portable endoscope with a built-in 4" high-definition screen – the game-changer for equine veterinarians on the race track! Designed for efficiency, this compact solution offers unparalleled portability, allowing veterinarians to perform on-the-spot examinations without the need for a dedicated facility. With a small diameter and advanced imaging capabilities, our endoscope ensures precise, detailed visuals of delicate areas, enabling thorough inspections without compromising the horse's comfort. What sets us apart is the ability to minimize or eliminate sedation, optimizing safety and expediting the diagnostic process. The real-time imaging on the 4" screen allows for immediate assessment, facilitating quick decision-making and efficient communication with horse owners and trainers. Cumbersome equipment is a dinosaur of the race track. Our battery powered endoscope is a cost-effective, time-saving tool that enhances your diagnostic capabilities and elevates the standard of care for racehorses. Your time matters, and so does theirs.

MODEL	DIA.	CHANNEL SIZE	LENGTH	DEFLECTION	FIELD OF VIEW
PB-953L	5.8mm	2.8mm	1500mm	180° U, 130° D, 120° L&R	120°

## See It All In High-Definition...

Our game-changing scope comes equipped with built-in video and image capture capabilities, ensuring that you can document and review every detail of your examinations effortlessly and directly from the buttons on the scope. The convenience continues with a compact 4" high-definition screen, seamlessly attached to the scope for immediate, crystal-clear visuals during on-the-go inspections. For those who prefer a larger display, we offer the flexibility of connecting to a 10" tablet-sized monitor, enhancing your viewing experience without compromising portability. Whether you're performing detailed examinations of delicate areas or capturing crucial footage for reference, our endoscope ensures that you have the advanced imaging tools you need, right at your fingertips. Elevate your diagnostic capabilities and streamline your workflow with our portable endoscope – where innovation meets efficiency for the modern equine veterinarian.



4" Screen Offers 720p High-Definition\*  
 10" Screen Offers Up To 1080p High-Defintion

A 3-Chip Camera Head Is Also Available With The  
 10" Display For Use With Dental Rigid Scopes.

