



TYPHOON

AUTOMATIC ENDOSCOPE REPROCESSOR



The All-New **TYPHOON AER** From **Rutledge Medical** Is The Latest Step In *Bio-Security Automation* For Veterinarians. Available In A Single Scope Or Dual Scope Basin, The Unit Can Complete An Entire Reprocessing Cycle In **15 Minutes**, Including The Leak Test, Rinse, Air Purge And Alcohol Purge Through All The Channels. Connections Are Available For All **Makes And Models** Of Endoscopes.

Product Model	AER-1
Machine size	550×660×940mm
Internal size	465×390×142mm
Cleaning quantity	1-Group (20 Min/Times)
Materials	ABS+Cold-Rolled Plates, PMMA, 304 Stainless Steel
Power Supply	220V 50HZ
Power	600W
Weight	75KG



✓ **Automatic Leak Testing**
Not just before the cycle, but during the cycle as well.
With safety features to emergency drain and terminate the cycle.

✓ **Sterile Water Rinsing**
Sterilized through three multi-micron filters as well as the germ eliminating power of the sun.

✓ **Multi-chemical Compatability**
Use various brands containing OPA, glutaraldehyde, hydrogen peroxide, paracetic acid and aldahol.

✓ **Hands-free Insert and Retrieve**
Load and unload your scope without contaminating the exterior of the unit, just use your foot.

✓ **Documentation**
An optional scope tracking system is available and the unit will provide a print out for every cycle with important data and verification.



AER-1
Automated Endoscope
Reprocessing System