



No Repairs, Just Results.



*EndoFresh® Veterinary portable endoscopy
system with disposable endoscopes*

Why Going Disposable?



Effective and maintenance-free

Every scope comes individually packaged and fully tested, ensuring consistent, top-tier performance. Everything will be predictable — no more hidden repair fees or last-minute emergency fixes.

Care on-the-go

Our portable video processor is designed for convenience. It packs into a compact suitcase, holding both the full system and your most-used endoscope, so you can easily take it wherever your care is needed.



Start, or resume your endoscopy practice today

Get your own practice up and running faster, at a fraction of the upfront cost. Eliminate the need for additional referrals and start treating your furry patients sooner, so they are back on their paws in no time!

Ask for tailored solutions

Every practice is different. That's why we offer customizable options to suit your unique needs – whether you're treating exotic species or diving into research. We'll work with you to create a solution that fits your practice perfectly.



Video Processor

Model	VCP-33G
Display	19 Inches
Signal	HDMI
Image storage	1TB
Report system	Yes
Integrated pump	Yes
Weight	15.2 KG
Dimension	440×375×149.5 mm

VCE Series Disposable Veterinary Endoscopes

Model	VCE-06S	VCE-15M	VCE-33L
Field of view (°)	120	140	140
Distal end dia. (mm)	5.4	8	13
Working length (mm)	600	1500	3300
Working channel (mm)	2.2	2.8	3.8
Angulation (°)	U/D: 210/210	U/D: 210/180 R/L: 160/160	U/D: 210/180 R/L: 160/160



Daichuan Medical (Shenzhen) Co., Ltd
+86 (0) 755-86576106
info@endofresh.com
Shenzhen Biomedical innovation Industrial Park
Pingshan, Shenzhen, China



Exclusively distributed in the USA by Rutledge Medical
+1 (888) ENDOVET
sales@endoscopes.vet
2828 E 51st Street, Suite 302, Tulsa, Oklahoma 74105, USA