



TYPHOON

AUTOMATIC ENDOSCOPE REPROCESSOR



The All-New **TYPHOON AER** From **Rutledge Medical** Is The Latest Step In *Bio-Security Automation* For Veterinarians. Available In A Single Scope Or Dual Scope Basin, The Unit Can Complete An Entire Reprocessing Cycle In **15 Minutes**, Including The Leak Test, Rinse, Air Purge And Alcohol Purge Through All The Channels. Connections Are Available For All **Makes And Models** Of Endoscopes.

Product Model	AER-1
Machine size	550×660×940mm
Internal size	465×390×142mm
Cleaning quantity	1-Group (20 Min/Times)
Materials	ABS+Cold-Rolled Plates, PMMA, 304 Stainless Steel
Power Supply	220V 50HZ
Power	600W
Weight	75KG



TYPHOON
UNO



Automatic Leak Testing

Not just before the cycle, but during the cycle as well.
With safety features to emergency drain and terminate the cycle.



Sterile Water Rinsing

Sterilized through three multi-micron filters as well as the germ eliminating power of the sun.



Multi-chemical Compatability

Use various brands containing OPA, glutaraldehyde, hydrogen peroxide, paracetic acid and aldahol.



Hands-free Insert and Retrieve

Load and unload your scope without contaminating the exterior of the unit, just use your foot.



Documentation

An optional scope tracking system is available and the unit will provide a print out for every cycle with important data and verification.



AER-1
Automated Endoscope
Reprocessing System