

Ultimate Tower Endoscopy system

VET-OR1800

Processor

Y Nr	Model Nr.	Description
Y100-263	VET-OR1800P	4K touch screen video processor w. 4K rigid camera
Y100-264	VET-OR1800	4K touch screen video processor w. FHD rigid camera
Y100-265	VET-OR1800E	4K touch screen video processor w/o rigid camera

Spare Parts

Y Nr	Model Nr.	Description
Y100-266	VET-OR1800 4K Camera Head Kit (Zoom Coupler)	4K camera head w. zoom coupler (14-32mm), light guide bundle, leakage tester&water resistant cap
Y100-0259	VET-OR1800 FHD Camera Head Kit (Zoom Coupler)	VET-CH160 w. zoom coupler (14-32mm), light guide bundle, leakage tester&water resistant cap

Trolley

Y Nr	Model Nr.	Description
TBD	AET-L1	4 shelves Trolley w. Monitor Arm, Holder for flexible endoscope and Rigid camera

What is Tower system?

- Image processor with light source, air and water pump,
- Work with flexible endoscope((VET-OR1800, VET-OR1200RH,VET-OR1200HD)
- Work with rigid endoscope (VET-OR1800,VET-OR1200RH)
- Need additional monitor and trolley
- Suitable for operating room use






3 LED

Improved clinical performance with 3 LED light source and 4K image processor



3 CBI

Target suspicious lesions and make accurate biopsy with 3 CBI



DICOM

Streamlined patient management with built-in DICOM



Flexible Endoscope for
VET-OR1800/ OR-160/ OR-100

Y Nr.	Type	Model	Working length (mm)	Diameter of distal end (mm)	Diameter of instrument channel (mm)	Range of bending
Y100-0235	Reusable	VET-2812R	1000	2.8	1.2	U180°D130°
Y100-0236	Reusable	VET-6028R	1500	6.0	2.8	U210°D90°L100°R100°
TBD	Reusable	VET-7828R	1500	7.8	2.8	U210°D90°L100°R100°
Y100-0217	Reusable	VET-FHD15	1500	9.8	2.8	U210°D90°L100°R100°
Y100-0218	Reusable	VET-FHD33	3300	9.8	2.8	U210°D90°L100°R100°
Y100-0219	Reusable	VET-FHD35	3500	12.8	3.7	U180°D180°L160°R160°
Y100-262	Adaptor, connect DEV Series to Veterinary Endoscopy System					
TBD	Disposable	DEV-3112	600	3.1	1.2	U210°D210°
TBD	Disposable	DEV-3822	600	3.8	2.2	U210°D210°
TBD	Disposable	DEV-4932	600	4.9	3.2	U210°D210°
Y100-260	Disposable	DEV-4922H	600	4.9	2.2	U210°D210°
Y100-261	Disposable	DEV-5832H	600	5.8	3.2	U210°D210°



Flexible Endoscope for
VET-OR1800/ OR-160/ OR-100

Canino

VET-6028R / VET-7828R / VET-FHD15

VET-2812R/DEV Series VET-2812R

Felino

VET-6028R / VET-7828R

VET-2812R/DEV Series VET-2812R

Equino

VET-FHD15 VET-FHD33 / VET-FHD35

VET-FHD15 VET-FHD33 / VET-FHD35

Delfines

VET-FHD33 / VET-FHD35

VET-FHD33 / VET-FHD35

